

NANOMATERIAL ELECTRICAL AND MECHANICAL PROPERTIES MEASURED IN A TRANSMISSION ELECTRON MICROSCOPE

Dmitri Golberg

Nanoscale Materials Centre, National Institute for Materials Science (NIMS),
Namiki 1-1, Tsukuba, Ibaraki, 3050044, Japan

World Premier International Centre for Materials Nanoarchitectonics (MANA), National Institute for Materials
Science (NIMS), Namiki 1-1, Tsukuba, Ibaraki, 3050044, Japan

E-mail: GOLBERG.Dmitri@nims.go.jp

I will demonstrate the significant progress in electromechanical property measurements of diverse low-dimensional individual nanostructures, e.g. nanotubes, nanowires, nanobelts and nanocones, carried out inside a transmission electron microscope (TEM). The measurements were performed by means of two high-resolution field-emission instruments operating at 300 kV, namely JEOL-3000F and JEM-3100FEF, and on free-standing individual objects using specialized piezo-driven scanning tunnelling microscope (STM)-TEM and atomic force microscope (AFM)-TEM holders (with conductive or non-conductive cantilevers).

A sharp etched Au tip, acting as the STM tip, or a Si cantilever, as the AFM probe, was assembled within the holders. Individual low-dimensional nanostructures were mounted on thin metallic (Au, Al) wires (\varnothing 250 μm), Fig. 1 and Fig. 2, which then were driven in three dimensions (with a precision better than 1 nm) inside the pole-piece of the microscopes using a piezo-motor until a tight physical contact between a structure and the probe was achieved.

The successful experimental runs to be discussed include multi-walled BN (Figs. 1 and 2) [1-4] and N-doped C nanotubes (NTs) [5,6], ZnO nanorods and nanowires [7], double-walled CNTs, multi-walled CNTs filled with Cu [8], Au, CuI or ZnS, inorganic nanothermometers, e.g. liquid Ga-filled MgO NTs and In-filled SiO₂ NTs, CdS nanobelts, and BN nano- and microcones.

The current-voltage (STM-TEM, Fig. 1) and/or force-displacement (AFM-TEM, Fig. 2) curves were recorded under a full control of the nanomaterial morphological, crystallographic and chemical changes/modifications during all stages of probing/manipulation. This allowed us to clearly interpret many new intriguing physical, chemical and electromechanical phenomena peculiar to novel advanced nanostructures. These include, but are not limited to deformation-driven electrical transport and piezoelectricity in individual Boron Nitride NTs, time-resolved electrical transport in N-doped CNTs, rheostat-like behaviour, femtogram (Cu) and attogram (CuI) mass transports, nanopipetting in Cu- and Cu-halide-filled CNTs, and strength and elasticity modifications of nanoscale composite materials made of filled C nanotubes.

References

1. D. Golberg et al., *Adv. Mater.* 19, 2413 (2007).
2. D. Golberg et al., *Nano Lett.* 7, 2146 (2007).
3. D. Golberg et al., *Acta Mater.* 55, 1293 (2007).
4. X.D. Bai, D. Golberg et al., *Nano Lett.* 7, 632 (2007).
5. P.M.F.J. Costa, D. Golberg et al., *Appl. Phys. A* 90, 225 (2008).
6. P.M.F.J. Costa, D. Golberg et al., *Appl. Phys. Lett.* 91, 2231081 (2007).
7. P.M.F.J. Costa, D. Golberg et al., *J. Mater. Sci.* 43, 1460 (2008).
8. D. Golberg et al., *Adv. Mater.* 19, 1937 (2007).
9. The author is indebted to Drs. P.M.F.J. Costa, X.D. Bai, M. Mitome, Y. Bando, and O. Lourie for their experimental and technical contributions.




Figure 1. (a) Consecutive TEM images of an individual multi-walled BN nanotube bent inside HRTEM; and (b) the corresponding I - V curves recorded at different stages of the experiment (a_1 - a_4). A notable current may pass through the nanotube in a bent state, while it disappears if the nanotube is straight; from Ref. [4].




Figure 2. Force-displacement curves directly recorded on a relatively thick (a) and a thin (b) individual multi-walled BN nanotube using an AFM-TEM setup. The insets display the nanotube appearance at different stages of the experiments: from Ref. [2].