



universität wien

Peter A. Lieberzeit
 Department of Analytical Chemistry
 and Food Chemistry
 Peter.Lieberzeit@univie.ac.at


Nanostructured Materials with Biomimetic Recognition Abilities for Chemical Sensing

Trends in Nanotechnology Conference, Tenerife, November 22nd, 2011




universität wien

What is a Sensor?




Data acquisition Measuring electronics Transducer Analyte Receptor

P. A. Lieberzeit, TNTConference, Tenerife, Nov. 22nd, 2011


universität wien

Manufacturing QCM




Basis: Quartz blanks (AT-cut), 168 μm thick

- ✓ Cleaning
- ✓ Screen-printing of the desired electrode pattern
- ✓ Burning of the electrode structures at 400°C

P. A. Lieberzeit, TNTConference, Tenerife, Nov. 22nd, 2011


universität wien

QCM - The Sauerbrey Equation



$f_0 = v/2d$

$$\frac{\Delta f}{f_0} = \frac{\Delta d}{d} = \frac{\Delta m}{r A d}$$

$$\Delta f = c \cdot f_0^2 \cdot \Delta m$$

P. A. Lieberzeit, TNTConference, Tenerife, Nov. 22nd, 2011


universität wien


Motivation

- Natural systems are outstanding in terms of **recognition and regeneration**.
- **But:** **limited ruggedness**, sometimes high cost.
- **Therefore:** Implement **biological functionality** into artificial materials.

P. A. Lieberzeit, TNTConference, Tenerife, Nov. 22nd, 2011








universität wien







Recognition Strategies




Outer plane – „concave“ Incorporation – „convex“

P. A. Lieberzeit, TNTConference, Tenerife, Nov. 22nd, 2011









Bioreactor Measurements

Bioreactor


- Cultivation of microorganisms under optimal conditions
- important parameters: the composition of the nutrient, oxygen supply, sterility, pH and temperature-controlled and controllable
- Sensors measure these parameters directly in the exhaust air or in the culture medium
- Process flow assessed









Lectins

- Glycoprotein
- Specific binding to oligosaccharide
 - Mannose
 - ConA, LCH, GNA
 - Galactose/N-acetylgalactosamine
 - RCA, FNA, AIL, VVL
 - N-acetylglucosamine
 - WGA, SNA, MAL
 - Fucose
 - UEA, AAL
- Play important role in cells and viruses recognition.

<http://www.ebi.ac.uk/pdbe-srv/ViewEntry/Zuvc/quaternary.html>





 **Acknowledgements**

S. Bajwa R. Samardzic T. Wangchareansak
G. Mustafa

 Austrian Science Fund (FWF) 

European Commission: EU-FP6 and FP7   

Royal Golden Jubilee Grant
PhD Program 

P. A. Lelercq, TINT-Conference, Tenerife, Nov. 22nd 2014

